

View Online at <https://aerobasegroup.com/nsn/6550-01-292-3743>

**Component Quantity:**

6

**Special Features:**

0consists of iron and transferrin in human serum; set contains 2 vials each of low, medium, and high level controls with lot specific label values for iron and tbc concentrations given in micrograms per deciliter; the controls should read within the following ranges of label values: plus/minus 15%, plus/minus 10%, plus/minus 10% respectively; sufficient for 20 runs of the low, medium, and high levels of iron concentration and tbc lot specific label values; 0.1% sodium aide as preservative; note: g

**Shelf Life:**

N/a

**Unit Of Measure:**

--

**Demilitarization:**

No

**Fiig:**

A23800