

View Online at <https://aerobasegroup.com/nsn/6550-01-292-3743>

Component Quantity:

6

Special Features:

0consists of iron and transferrin in human serum; set contains 2 vials each of low, medium, and high level controls with lot specific label values for iron and tbc concentrations given in micrograms per deciliter; the controls should read within the following ranges of label values: plus/minus 15%, plus/minus 10%, plus/minus 10% respectively; sufficient for 20 runs of the low, medium, and high levels of iron concentration and tbc lot specific label values; 0.1% sodium aide as preservative; note: g

Shelf Life:

N/a

Unit Of Measure:

--

Demilitarization:

No

Fiig:

A23800