Blood Gas And Ph Control - Page 1 of 1



View Online at https://aerobasegroup.com/nsn/6550-01-329-9524

Physical Form:

Liquid

Usage Design:

For monitoring the measurement of ph, pco2 and po2 in blood gas instrumentation and for in vitro diagnostic use

Primary Container Type:

Ampul

Primary Container Content Quantity:

1.700 milliliters

Total Primary Container Quantity:

50

Special Features:

Level 2 - normal (yellow ampul), buffered bicarbonate in equilibrium with pre determine levels of oxygen, carbon dioxide and nitrogen

Shelf Life:

N/a

Unit Of Measure:

--

Demilitarization:

No

Fiig:

A10000