

View Online at <https://aerobasegroup.com/nsn/6550-01-487-2904>

Component Breakdown Within Each Test Kit:

96 cytomegalovirus antigen coated microwell strips; 1 50 ml sample diluent bottle; 11 ml enzyme conjugate; 0.3ml each vial negative calibrator reference standard set; 11 ml peroxide substrate solution a; 11 ml tmb chromogen solution b; 1 vial 50 ml washing buffer (20x); 1 well holder

Shelf Life:

N/a

Unit Of Measure:

--

Demilitarization:

No

Fig:

A510v0