

View Online at <https://aerobasegroup.com/nsn/6550-01-487-2909>

Component Breakdown Within Each Test Kit:

96 rubella antigen coated microwell strips; 1-50 ml sample diluent bottle; 11 ml anti-human igg conjugated to horseradish peroxidase; 2 ml each vial negative reference standard set; 11 ml peroxide substrate solution a; 11 ml tmb chromogen solution b; 1-50 ml vial washing buffer; one each well holder

Shelf Life:

N/a

Unit Of Measure:

--

Demilitarization:

No

Fig:

A510v0