

View Online at <https://aerobasegroup.com/nsn/6550-01-493-2606>

**End Application:**

Reagent cartridges for performing cardiac testing profiles.

**Product Name:**

Cardiac panel triag 25s

**Special Features:**

The disposable cardiac panel measures the concentration of creatine kinase mb (ck-mb), myoglobin and troponin i; the panel is stable while in the sealed pouch when stored 2 to 8 degrees c (35 to 46 degrees f); once removed from refrigeration the panel is stable up to 14 days; the manufacturer is biosite diagnostic part number 97000

**Shelf Life:**

N/a

**Unit Of Measure:**

--

**Demilitarization:**

No