

View Online at <https://aerobasegroup.com/nsn/6550-01-500-2383>

**End Application:**

For use in performance testing of culture media, stains, biochemical id kits, maintenance of stock cultures and in the evaluation of bacteriological procedures

**Product Name:**

Culture control anaerobic panel microscan multi-loop set 20s

**Special Features:**

Microscan multi-loop set; to grow anaerobic culture control organisms for quality control of microscan instruments; contains stabilized, preserved viable microorganisms; each loop may be dissolved in re-hydration media or plated/streaked directly onto appropriate media; includes: (5) bacteroides fragilis atcc 25285, (5) clostridium perfringens atcc 13124, (5) clostridium sordellii atcc 9714, and (5) peptostreptococcus magnus atcc 29328; should be stored upon arrival 2 to 8 degrees c (or frozen for campylobacter sp.) each loop is individually packaged in a foil pouch with a peel-off label

**Shelf Life:**

N/a

**Unit Of Measure:**

--

**Demilitarization:**

No