

View Online at <https://aerobasegroup.com/nsn/6550-01-573-9356>

Storage Temperature Range:

+2.0/+8.0 degrees celsius

Product Name:

Snap 4dx test 30s

Special Features:

Screening test for canine heartworm, e. Canis, lyme disease, and a. Phagocytophilum; kit contains 8ml bottle anit-hw/ap/ly/ec: hrpo conjugate, reagent rack, transfer pipets, sample tubes, and snap devices; store at 2 to 8 degrees c until expiration date; after snap devices and reagents are removed from 2 to 8 degrees c for more than 24 hours, expiration date is 90 days or the printed expiration date, whichever comes first; package of 30

Shelf Life:

N/a

Unit Of Measure:

--

Demilitarization:

No