NSN 6550-01-644-2572

Dna Analysis Test Kit - Page 1 of 1

A510v0



View Online at https://aerobasegroup.com/nsn/6550-01-644-2572

Volume at https://doi.obdoog.cap.com/mon/cocc c.f. c.f.f. 2012
Fotal Primary Container Quantity:
100
Product Name:
ndicating fta microcard, 100s
Test Media:
Blood
Component Breakdown Within Each Test Kit:
Wb 1000003, desiccant; wb 100036, multi-barrier pouch, small; wb 120204, fta purification reagent; wb 100030, fta gene card trays; w
100032, sterile foam tipped applicator swabs
Special Features:
Eta indicating micro cards have a color indicator that changes from pink to white when sample is applied to prevent false positive sampling;
recommended for use with clear samples such as buccal or cultured cells; fta micro cards allow for samples of 125?L whole blood or 25?L
plant homogenate within the sample area; capture nucleic acids in one easy step; simply apply your sample to the fta card; nucleic acids
remain immobilized and are stabilized for transport, immediate processing or long-term room temperature storage; applications: genetic
dentification; animal breeding studies; molecular biology; transgenic identification, plasmid screening, drug discovery, genomics, whole
genome amplification; genomic dna stored on fta cards may be tored at room temperature for more than 17.5 years; pack of 100
Shelf Life:
N/a
Unit Of Measure:
-
Demilitarization:
No
Filg: