

View Online at <https://aerobasegroup.com/nsn/6550-01-650-4412>

**Product Name:**

Snap 4dx plus test, 30s

**Special Features:**

Snap 4dx plus test; canine screening for 6 vector-borne diseases (e.Canis, e. Ewingii ab, lyme disease ab, a. Phagocytophilum/a. Platys ab, and heartworm ag) in canine serum, plasma or whole blood; kit c/o reagent, wash solution 0.4 ml, substrate solution 0.6 ml and 30 snap devices; store at 2 to 8 degrees c; the snap 4dx plus test requires 3 drops of sample and 4 drops of conjugate; snap devices and reagents can be stored at room temperature 18?25?C for 90 days or until the expiration date, whichever occurs first

**Shelf Life:**

N/a

**Unit Of Measure:**

--

**Demilitarization:**

No