

View Online at <https://aerobasegroup.com/nsn/6630-01-046-7124>

**General Description:**

Operating modes purge, desorb, trap bake; trap temp control direct set proportional controller, ambient to 250 degrees c; purge timer 0-30 min. Variable; purge valve solenoid valve with automatic shutoff; manual purge/desorb valve; 19 in. H, 21 in. W, 15 in. D; 115 vac, 60 hz, 2 amp; 5 ml sampler, 5 ml sample syringe, 10 ul calibration syringe, 5 spare septa, spare o-rings and ferrules, blank spare trap column, tube/needle sample transfer system, fuses, and manual; concentrates organic compounds in fresh and waste water

**Shelf Life:**

N/a

**Unit Of Measure:**

--

**Demilitarization:**

No