

View Online at <https://aerobasegroup.com/nsn/6630-01-268-8452>

Special Features:

Test kit bath oxidation 115/230 vac, 50/60 hz, interior illuminated through frosted pyrex safety shield, oxygen flowmeters controlled by needle valve at each tube, temperature controlled at 95 +/- .2 degrees c. By solid state controller and thermister, insulated pyrex flask, agitated by a turbine stirrer, accomadates 8 tubes

Shelf Life:

N/a

Unit Of Measure:

--

Demilitarization:

No

Fig:

A23800