

View Online at <https://aerobasegroup.com/nsn/6630-01-529-6134>

End Application:

Stat profile critical care xpress - ccx 2, 2+, or ccx 1+ configured as ccx 2/2+

Test Type For Which Designed:

Blood gasses - ph, pco2, po2; so2 percent; electrolytes - na+, k+, ci-, ca++, mg++; chemistry - gluc, bun, creat, lact; hematology - hct, hb

Test Sample/specimen Type:

Whole blood

Special Features:

Part name assigned by controlling agency - ccx 2 calibration cartridge, with creatinine

Shelf Life:

N/a

Unit Of Measure:

--

Demilitarization:

No

Fiig:

A503u0