NSN 6630-01-537-9962

Water Testing Swab - Page 1 of 1



View Online at https://aerobasegroup.com/nsn/6630-01-537-9962

Product Name:

Swab water atp watergiene 100s

Special Features:

Watergiene detects total atp in water, which is a critical indicator of total biologic activity resulting from microorganisms, as well as plant, animal or food material; the watergiene can be used with the novalum, nsn: 6630-01-536-1456, or with the lumgiene, nsn: 6640-01-530-1441; item needs to be refrigerated; detects atp (adenosine triphosphate) in water in seconds; a simple twist motion activates reagents before the swab is inserted into and read on a luminometer; box/100

| 0 | L - | . I£ | ifc | |
|---|-----|------|-----|--|
| | | | | |

N/a

Unit Of Measure:

--

Demilitarization:

No