

View Online at <https://aerobasegroup.com/nsn/6640-01-340-9658>

**Component Quantity:**

2

**Nonsupply Items And Quantities:**

15 ml lactic acid in stopper btl w/cap 1; 30 ml titration reagent in stopper btl w/cap 1

**Special Features:**

Direct reading titration method; 0-200 ppm cl<sub>2</sub> in 4 ppm increments; performs starch iodide-thiosulfate titrate analysis; epa accepted; sufficient qty to perform 50 test at 200 ppm

**Shelf Life:**

N/a

**Unit Of Measure:**

--

**Demilitarization:**

No

**Fiig:**

A23800