

View Online at <https://aerobasegroup.com/nsn/6640-01-340-9658>

Component Quantity:

2

Nonsupply Items And Quantities:

15 ml lactic acid in stopper btl w/cap 1; 30 ml titration reagent in stopper btl w/cap 1

Special Features:

Direct reading titration method; 0-200 ppm cl₂ in 4 ppm increments; performs starch iodide-thiosulfate titrate analysis; epa accepted; sufficient qty to perform 50 test at 200 ppm

Shelf Life:

N/a

Unit Of Measure:

--

Demilitarization:

No

Fiig:

A23800