

View Online at <https://aerobasegroup.com/nsn/6640-01-652-9923>

**Special Features:**

Hydrogen generator for use with gas chromatograph mass spectrometer (gcms) products. Produces hydrogen gas to operate gcms. Kit includes power cord, user manual, silica gel, deionizer bags, desiccant cap tool, 1/8 in refrigeration-grade copper line with 1/4 in and 1/8 in swagelok fittings for connection to the griffin gc/ms, and a fitting that connects the generator to a 0.375 in (3/8 in) tygon tube, which is inserted into a deionized water container or water supply. Ref: 664001d1. H9d/socom logistics help desk: cml (859)293-3114; dsn: 745-3114; toll free: 888-763-7259; email: [logisticshelpdesk@sofsa.Mil](mailto:logisticshelpdesk@sofsa.Mil).

**Shelf Life:**

N/a

**Unit Of Measure:**

--

**Demilitarization:**

No

**Fiig:**

T396-g