

View Online at <https://aerobasegroup.com/nsn/6665-01-467-3718>

Component Quantity:

17

Product Name:

Water testing kit bacteriological 17components mel/mf total coli

Special Features:

Mel/mf total coliform lab includes usepa accepted m-endo broth pour rite ampules for detecting total coliform contaminations within 24 hours.

The endo broth can be used to test waste, surface, ground and recreational waters and saltwater, provided the turbidity is low enough.

Components: step by step procedures manual, 200 m-endo broth pour rite ampules, portable incubator, mf petri dish rack for 42 petri dishes, 200 whirl-pak bags w/dechlorinating agent, field filtration apparatus w/funnels and filters, 200 petri dishes w/pads, and alcohol burner, armoured alcohol thermometer, 5 germicidal cloths

Shelf Life:

N/a

Unit Of Measure:

--

Demilitarization:

No

Fiig:

A23800