

View Online at <https://aerobasegroup.com/nsn/6840-00-926-9117>

**Physical Form:**

Liquid

**Product Type:**

Glutaraldehyde

**Composition Percentage By Weight:**

2.0 glutaraldehyde

**Test Organism Type:**

Bacillus subtilis and escherichia coli and herpes simplex influenza vaccinia and klebsiella pneumoniae and pseudomonas aeruginosa, atcc 10145 and staphylococcus aureus, atcc 6538 and streptococcus pyogenes and trichophyton interdigitale, atcc 640

**Dilution By Volume For Testing:**

Undiluted

**Test Time Limit:**

180.0 minutes and 1.0 minutes and 1.0 minutes and 1.0 minutes and 1.0 minutes and 1.0 minutes and 1.0 minutes and 0.5 minutes

**Test Temp:**

20.0 degrees celsius

**Product Name:**

Disinfectant general purpose glutaraldehyde liquid 1 gallon

**Special Features:**

When activated by an accompanying vial w/dye produces a 7.0 to 8.9 ph solution as a rapid-acting (not more than 2 hrs) sporicidal, tuberculocidal, virucidal, bactericidal, fungicidal, and pseudomonacidal agent; corrosion inhibiting

**Supplementary Features:**

Use within 14 days after activation; (3.78 liters); suitable for disinfecting and cold sterilization of surgical instruments

**Shelf Life:**

N/a

**Unit Of Measure:**

--

**Demilitarization:**

No

**Fig:**

T090-b